



## Medical Ozone Generator Accessories

[www.aquapureozone.com](http://www.aquapureozone.com)

# Oxygen Regulator



AOT-OR-870



AOT-OR-540



AOT-OR-BN1

AOT-OR-BN2



Oxygen flow setting: 1/32, 1/16, 1/8, 1/4, 1/2, 3/4, 1, 1.5, 2, 3, 4 L/min  
Tube connector: Barb type

# Ozone Destructor

## Product Description :

This catalyst is used to effectively destroy ozone emitted from various off-gas emissions, converting ozone to oxygen since it is a true catalyst it is not consumed during the destruction of the ozone, so its lifetime is extensive. Inlet ozone concentration has no bearing on the design of the destruct system with our catalyst.

**Connector** : PVDF material

**Particle size** : 4.8 x 2.4 mm

**Density** : 850 kg/m<sup>3</sup>

**O3 input** : Normal connector/ luer connector

**Storage condition** : Cool and dry area

## AOT-OD-20

Product size: 215x18x18 MM  
Weight: 0.2kg

## AOT-OD-10

Product size: 115x18x18 MM  
Weight: 0.1kg



# Accessories

## Syringe/Luer lock Set



### Specification:

- ◆ 30ml
  - ◆ 60ml
  - ◆ 100ml
  - ◆ 200ml
- With luer lock connector



### Specification:


With male and female  
Material: PP  
Recommend tube size: 4-5mm(inner diameter)  
Size: 3/16" (4.8mm outer diameter )

# Diffuser




## Specification:

- Inlet connector: OD 5mm
- Low cost
- Easy to set up
- For water treatment
- Need air pump



# Accessories

## Check valve



### Specification:

Size: 36x19mm  
Connector: OD 5.8mm  
For air flow: <1.4L/min  
Material: Viton (PVDF available)  
Weight: 0.002kg

## Silicone Tube



### Specification:

FDA standard  
Size: 7x4mm or 8x4 mm  
Thickness: 1.5mm  
Material: Silicone (FDA)  
Color: transparent  
Ozone resistance

## Luer filter



### Specification:

OD: 20MM  
Pore size: 0.22  
Material: ABS  
Outlet: Male Luer Lock with Thread  
Inlet : Female Luer Lock with Thread  
For both water and air



## Connector



2 ways  
OD4.8mm(3/16")  
PP material




3 ways  
OD4.8mm(3/16")  
PP/PVDF material



4 ways  
OD4.8mm(3/16")  
PP material



5 ways  
OD4.8mm(3/16")  
PP material



6 ways & 8 ways  
Main port OD:8mm  
Shunt port OD: 4mm  
PP material



Oxygen Inlet port



## Specification:

Size: 40.5 x 11 mm  
Diameter: OD5mm, ID3.4mm  
Material: 316 stainless steel  
Weight: 8g

Ozone outlet port



## Luer lock specification:

Material : 316 Stainless Steel  
Inner Diameter : 5.0 mm (with thread) and 2.4 mm  
Outer Diameter : 7.8 mm (with thread) and 5.0 mm  
Length : 29.4 mm  
Weight : 5g




# Accessories

## Stopcock



### Specification:

Material : PC  
With male and female luer connector



# Accessories

## Ozone bag/Ozone glass cup




### Specification:

Material : Teflon  
Connector: PC material  
Volume: 500ml  
Working temperature: -240°C-+250°C  
Connection size: 3/16




### Specification:

Volume: 500ml  
With stopcock  
Straight pipe : 15cm  
L pipe: 4.5\*10.5cm



# Accessories

## Luer catheter



### Specification:

Material : PVC

OD: 4.7mm

ID: 3.4mm

14 french



# THANK YOU!



## **AQUAPURE Ozone Technology Co., Ltd**

R1102, University Town Business Park, Lishan Road, Taoyuan, Nanshan, 518055 Shenzhen, China

Tell : 0086-755-2660 5963

Mobile: 0086-13715076653

Email: [info@aquapureasia.com](mailto:info@aquapureasia.com)